

Product sheet

OS-RC-2 | 305086

Flask Coating
The flask is coated with a special coating to ensure optimal cell attachment and growth.

Freezing Procedure
The cells should be frozen in a controlled manner, using a freezing medium and a freezing container, and stored at -78°C.

Shipping Conditions
The cells should be shipped in a controlled manner, using a shipping container, and stored at -78°C.

Storage Conditions
The cells should be stored in a controlled manner, using a storage container, and stored at -150 °C for 196 days.

HLA / HLA

Sterility
The cells are sterile and free of mycoplasmas. The PCR results are negative for all HLA markers.