





WI 38 VA13 ~~WI-2RA~~ ~~WI-2RA~~ | 300421

### Thawing and Culturing Cells

1. ~~...~~
2. ~~...~~
3. ~~...~~
4. ~~...~~
5. ~~...~~
6. ~~...~~
7. ~~...~~
8. ~~...~~

### Incubation Atmosphere

37°C, 5% CO<sub>2</sub>, ~~...~~

### Flask Coating

~~...~~

### Freezing Procedure

~~...~~ -78°C

### Shipping Conditions

~~...~~ -78°C

### Storage Conditions

~~...~~ -150 °C 196 ~~...~~

~~...~~ / ~~...~~ / HLA

### Sterility

~~...~~ PCR ~~...~~

~~...~~