

Product sheet

KLN-205 | 400419

DESCRIPTION

Description KLN-205 is a murine B cell hybridoma cell line that secretes monoclonal antibodies against the mouse IgG2b antibody. KLN-205 is a murine B cell hybridoma cell line that secretes monoclonal antibodies against the mouse IgG2b antibody.

Organism Mouse

Tissue B cell

Disease B cell lymphoma

Synonyms KLN 205, KLN205

CHARACTERISTICS

Breed/Subspecies DBA/2

Growth properties Adherent

IDENTIFICATION

Citation KLN-205 (ATCC CRL-2739) Cytion 400419

Biosafety level 1

NCBI_TaxID 10090

CellosaurusAccession CVCL_3533

TESTING

Tumorigenic No, DBA/2, BDF1

CULTURE

Culture Medium EMEM (MEM Eagle), w: 2 mM L-Glutamine, w: 2.2 g/L NaHCO3, w: EBSS (ATCC CRL-2739 Cytion 820100a)

