

Product sheet

NRK-IBB-DiHcRed1 | 500671

Incubation Atmosphere 37°C, 5% CO₂, humidified

Flask Coating Cell culture grade, non-adhesive

Freezing Procedure Harvest cells into 10% FBS serum, centrifuge at 300g for 5 min, resuspend in freezing medium, aliquot into 1 ml vials, store at -78°C

Shipping Conditions Dry ice, -78°C

Storage Conditions -150 to -196 °C, 196 °C, 196 °C

HLA

Sterility PCR, 196 °C