

XXXXA9 | 305166

XXXXXXXX XXXX-XXXXXXXXXXXX

Antigen expression H-2k

Tumorigenic XX, XXXXXXX XXXXXXX

XXXXXX

Culture Medium DMEM, w: 4.5 g/L XXXXXXX, w: 4 mM L-XXXXXXX, w: 3.7 g/L NaHCO₃, w: 1.0 mM XXXX XXXXXXX (XXXX XXXXXXX Cytion 820300a)

Supplements XXXX XXXXX 10% FBS

Dissociation Reagent XXXXXXX

Subculturing XXX XX XXXXXXX XXXX XXXXXXX XXXXXXX XXXXX XXXX X-PBS XXX XXXX XXXXXXX XXXX XXXXXXX T25, XXXXX X-3-5 X' X-PBS, XXXXX XXXX 3 XXXX. XXXX XX XXXXXXX XXXXXXX, XXXX XX XXXXXXX XXXXXXX XXXXXXX XXX XXXXXXX XXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXX.

Fluid renewal 2 X 3 XXXXXXX XXXXXXX

Freeze medium XXXXXXX XXXXXXX XXXXXXX, XXX XXXXXXX XXXXXXX XXXXXXX XXXX (XXXX FBS) + 10% DMSO XXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX XXXXXXX, XXX C

