

RPMI 1640, w: 2.0 mM , w: 2.0 g/L NaHCO3 |
820700a

XXXXX RPMI 1640, XXXXXX XX XXXXXX RPMI, XXX XXXXX XXXXXX XXXX X-XXXXXXXX XXXXXX XXXXX XXXXXX XXXXXX XXXXXX XXXXXX XXXXXX XXXXXX XXXXXX XXXXXX

Media: RPMI 1640 Medium (Gibco) supplemented with 10% FBS (Gibco), 100 U/ml penicillin (Gibco), 100 U/ml streptomycin (Gibco), and 100 U/ml nystatin (Gibco).

RPMI 1640 (Gibco) supplemented with 10% Fetal Bovine Serum (FBS) (Gibco) and 100 units/ml penicillin, 100 units/ml streptomycin, and 100 units/ml nystatin. Cells were grown in 25 cm² flasks (Corning) at 37°C in 5% CO₂.

Media: RPMI 1640 Medium + 10% FBS + 5-10% CO₂, pH 7.2-7.4

- [illegible]

		(" / ')
XXXXXXXX XXXX	XXXXXX	10.0
	L-XXXXX-L-XXXXXX	434.40
	L-XXXXXX	200.0
	L-XXXXXX H2O	56.82
	XXXXX L-XXXXXX	20.0
	L-XXXXXX 2HCl	65.2
	XXXXX L-XXXXXX	20.0
	L-XXXXXX HCl H2O	20.27

Product sheet

RPMI 1640, w: 2.0 mM
820700a, w: 2.0 g/L NaHCO3 |

	L-Asparagine-L-Asparagine	20.0
	L-Glutamine	50.0
	L-Proline	50.0
	L-Leucine HCl	40.0
	L-Isoleucine	15.0
	L-Valine	15.0
	L-Methionine	20.0
	L-Phenylalanine	30.0
	L-Histidine	20.0
	L-Tyrosine	5.0
	L-Glycine 2Na 2H2O	28.83
	L-Serine	20.0
Vitamins	Thiamine p-Ascorbic Acid	1.0
	D-Biotin	0.2
	Inositol Nicotinic	3.0
	D-Calcium Pantoic	0.2
	Nicotinic Nicotinic	1.0
	Pyridoxal-Hydrochloride	35.0
	Pyridoxal-Hydrochloride	1.0
	Pyridoxal-Hydrochloride HCl	1.0
	Pyridoxal-Hydrochloride	0.2

Product sheet

RPMI 1640, w: 2.0 mM
820700a , w: 2.0 g/L NaHCO3 |

	Hydrochloric acid	1.0
	Vitamin B12	0.0
Calcium nitrate	Ca(NO ₃) ₂ · 24H ₂ O	100.0
	KCl	400.0
	MgSO ₄ · 7H ₂ O	100.0
	NaCl	6000.0
	NaHCO ₃	2000.0
	Na ₂ HPO ₄	800.00
DL-phenylalanine	D-phenylalanine	2000.0
	L-phenylalanine	1.00
	DL-phenylalanine	5.30