

BJ | 305222

Karyotype BJ

Culture Medium DMEM w 4.5g/L w 4mM L- w 3.7g/L NaHCO3 w 1.0mM Cytion 820300a

Supplements 10% FBS 20 ng/mL bFGF

Dissociation Reagent

Subculturing PBS T25 3-5ml T75 5-10ml PBS T25 1-2m

Freeze medium FBS 10 DMSO CM-1

Thawing and Culturing Cells

- 1.
2. -150 3
3. 37 40 60
4. 70
5. 8ml 15ml
6. 300 x g 3
7. 10ml 2 T25 1
- 8.

Incubation Atmosphere 37 5% CO2

BJ | 305222

Flask Coating

Freezing Procedure -78

Shipping Conditions -78

Storage Conditions -150 -196 80

HLA

Sterility PCR