

**Ba/F3 | 305224**

<b>Description</b>	BALB/c	B	BA/F3	BaF3					
	Baf3	IL-3		Ba/F3	F3	Ba/F3	F3		IL-3
	BA/F3						CML	BCR-ABL	
		BA/F3							

<b>Organism</b>	
<b>Tissue</b>	

<b>Synonyms</b>	BA/F3 BaF3 BAF3 Baf3
-----------------	----------------------

<b>Breed/Subspecies</b>	C3H
-------------------------	-----

<b>Morphology</b>	
<b>Cell type</b>	B

<b>Growth properties</b>	
--------------------------	--

<b>Citation</b>	Ba/F3	305224
-----------------	-------	--------

<b>Biosafety level</b>	1
------------------------	---

<b>NCBI_TaxID</b>	10090
-------------------	-------

<b>CellosaurusAccession</b>	CVCL_0161
-----------------------------	-----------

<b>Karyotype</b>	Ba/F3	33
------------------	-------	----

**Ba/F3 | 305224**

**Culture Medium** RPMI1640 w 2.0mM w 2.0g/L NaHCO3 Cytion article number 820700a

**Supplements** 5 FBS 10 ng/mL mouse IL-3

**Subculturing** <sup>6</sup> /ml 5×10<sup>5</sup> /ml 3×10<sup>5</sup> 1×10

**Freeze medium** FBS 10 DMSO

CM-1

**Thawing and Culturing Cells**

- 1.
2. -150 3
3. 37 40 60
4. 70
5. 8ml 15ml
6. 300 x g 3
7. 10ml 2 T25
- 8.

1

**Incubation Atmosphere** 37 5% CO<sub>2</sub>

**Flask Coating**

**Freezing Procedure** -78

**Ba/F3 | 305224**

**Shipping  
Conditions**

-78

**Storage  
Conditions**

-150 -196

80

### **HLA**

**Sterility**

PCR