

HEK293-CTLA4 | 305423

Culture Medium RPMI1640 w 2.0mM w 2.0g/L NaHCO3 Cytion article number 820700a

Supplements 10% FBS 1mM 10mM HEPES 1% NEAA 1 mg/mL Geneticin G418-Sulfat

Dissociation Reagent -EDTA

Subculturing PBS PBS_{CO2} 2 3 Trypsin/EDTA T25

Fluid renewal 2 3

Post-Thaw Recovery T25 1 2 1 3 24

Freeze medium FBS 10 DMSO CM-1

Thawing and Culturing Cells

- 1.
2. -150 3
3. 37 40 60
4. 70
5. 8ml 15ml
6. 300 x g 3
7. 10ml 2 T25 1
- 8.

Incubation Atmosphere 37 5% CO₂

HEK293-CTLA4 | 305423

Flask Coating

Freezing Procedure -78

Shipping Conditions -78

Storage Conditions -150 -196 80

HLA

Sterility PCR