

HSF SV40 | 305338

Description	HSF(SV40)	40 SV40	T	T-Ag	SV40 T-Ag
	HSF(SV40)			SV40 T-Ag	

Organism

Morphology

Cell type

Growth properties

Citation	HSF(SV40) (Cytion	305338)
-----------------	-------------------	---------

Biosafety level 1

NCBI_TaxID	9606
-------------------	------

GMO Status GMO-S1 HSF SV40 T-

Culture Medium	DMEM:Ham's F12 (1:1) w: 3.1 g/L	w: 2.5 mM L-	w: 15 mM HEPES	w: 0.5 mM	w: 1.2 g/L NaHCO ₃ (Cytion)
-----------------------	---------------------------------	--------------	----------------	-----------	--

Supplements 10% FBS 50µg/ml

Dissociation Reagent

HSF SV40 | 305338

Freeze medium

FBS 10 DMSO

CM-1

Thawing and Culturing Cells

- 1.
2. -150 3
3. 37 40 60
4. 70
5. 8ml 15ml
6. 300 x g 3
7. 10ml 2 T25
- 8.

Incubation Atmosphere

37 5% CO₂

Flask Coating

Freezing Procedure

-78

Shipping Conditions

-78

Storage Conditions

-150 -196 80

HLA

HSF SV40 | 305338

Sterility

PCR