

TM4 | 305170

Antigen expression H-Y , Musculus,

Tumorigenic

Culture Medium DMEM:Ham's F12 (1:1) w: 3.1 g/L w: 2.5 mM L- w: 15 mM HEPES w: 0.5 mM w: 1.2 g/L NaHCO3 (Cytion)

Supplements 2.5% FBS 5%

Dissociation Reagent

Subculturing PBS T25 3-5ml T75 5-10ml PBS T25 1-2ml

Fluid renewal 2 3

Freeze medium FBS 10 DMSO CM-1

Thawing and Culturing Cells

- 1.
2. -150 3
3. 37 40 60
4. 70
5. 8ml 15ml
6. 300 x g 3
7. 10ml 2 T25 1
- 8.

