

HROC113 | 300803

Description	2006	PD	CRC
--------------------	------	----	-----

Organism

Tissue	UICC IV
---------------	---------

Disease TNM T4N2M0R0 G3 Lk(n)+5 Σ Lk(n)45

Synonyms	HROC113P
-----------------	----------

Age	41
------------	----

Gender

Ethnicity

Morphology

Growth properties

Citation	HROC113 (Cytion	300803)
-----------------	-----------------	---------

Biosafety level 1

NCBI_TaxID	9606
-------------------	------

CellosaurusAccession CVCL_S852

Protein expression

HROC113 | 300803

Tumorigenic

Viruses SV40 JC/BK HBV HCV HIV

Ploidy status

MSI-status MSI-H

Mutational profile P53 APC K-RasG12VA N-Raswt H-Raswt PIK3CAwt B-Rafwt

Culture Medium DMEM:Ham's F12 (1:1) w: 3.1 g/L w: 2.5 mM L- w: 15 mM HEPES w: 0.5 mM w: 1.2 g/L NaHCO3 (Cytion)

Supplements 10% FBS

Dissociation Reagent

Doubling time 23

Subculturing PBS T25 3-5ml T75 5-10ml PBS T25 1-2m

Seeding density 2×10^4 /cm²

Fluid renewal 3 5

Post-Thaw Recovery

Freeze medium FBS 10 DMSO CM-1

HROC113 | 300803

**Thawing and
Culturing Cells**

- 1.
2. -150 3
3. 37 40 60
4. 70
5. 8ml 15ml
6. 300 x g 3
7. 10ml 2 T25 1
- 8.

Incubation Atmosphere 37 5% CO₂

Flask Coating

Freezing Procedure -78

Shipping Conditions

-78

Storage Conditions -150 -196 80

HLA

Sterility PCR