

hGF | 300703

Description									
hGF		*Porphyromonas gingivalis*							
Organism									
Tissue									
Applications									
Cell type									
Growth properties									
Citation									
hGF		Cytion	300703						
NCBI_TaxID		9606							
Culture Medium									
DMEM:Ham's F12 (1:1)		w: 3.1 g/L	w: 2.5 mM L-	w: 15 mM HEPES	w: 0.5 mM	w: 1.2 g/L NaH			
Supplements		10% FBS 10 ng/mL bFGF 10 microgram/L Insulin							
Dissociation Reagent									
Subculturing									
PBS		T25	3-5ml	T75	5-10ml	PBS	T25		

hGF | 300703

Freeze medium	90	FBS	10	DMSO	CM-1	Cytion
---------------	----	-----	----	------	------	--------

Thawing and Culturing Cells	1.					
	2.			-150	3	
	3.			37	40	60
	4.	70				
	5.		8ml	15ml		
	6.	300 x g	3			
	7.	10ml		2	T25	1
	8.					

Incubation Atmosphere	37	5% CO ₂
-----------------------	----	--------------------

Flask Coating

Freezing Procedure	-78
--------------------	-----

Shipping Conditions	-78
---------------------	-----

Storage Conditions	-150	-196	80
--------------------	------	------	----

Product sheet

hGF | 300703

Sterility

PCR