

BJ | 305222

**Description**

BJ  
BJ  
BJ, hTERT  
, BJ BJ

**Organism**

**Tissue**

**Synonyms** FF-WT-BJ, BJ1

**Age** 1

**Gender**

**Ethnicity**

**Morphology**

**Cell type**

**Growth properties**

**Citation** BJ( 305222)

**Biosafety level** 1

**NCBI\_TaxID** 9606

**CellosaurusAccession** CVCL\_3653

**BJ** | 305222

**Karyotype** BJ . .

**Culture Medium** DMEM, w: 4.5g/L , w: 4mM L- , w: 3.7g/L NaHCO3, w: 1.0mM ( 820300a)

**Supplements** 10% FBS, 20ng/mL bFGF

**Dissociation Reagent**

**Subculturing** PBS . T25 3~5ml, T75 5~10ml PBS . T25

**Freeze medium** (FBS ) + 10% DMSO , CM-1(Cytion 800100)

**Thawing and Culturing Cells**

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

**Incubation Atmosphere** 37°C, 5% CO2' .

BJ | 305222

---

**Flask Coating**

**Freezing Procedure**

-78°C

**Shipping Conditions**

-78°C

**Storage Conditions**

-150°C -196°C .80°C

/ / HLA

**Sterility**

PCR