

Ba/F3 | 305224

Description	BALB/c	-B	BA/F3						
	Baf3		IL-3	. Ba/F3	F3	Ba/F3	. F3		IL-3
	BA/F3					. ,	(CML)	BCR-ABL	
	BA/F3								

Organism	
Tissue	

Synonyms	BA/F3, BaF3, BAF3, Baf3
-----------------	-------------------------

Breed/Subspecies	C3H
-------------------------	-----

Morphology	
Cell type	Pro-B

Growth properties	
--------------------------	--

Citation	Ba/F3(305224)
-----------------	----------------

Biosafety level	1
------------------------	---

NCBI_TaxID	10090
-------------------	-------

CellosaurusAccession	CVCL_0161
-----------------------------	-----------

Karyotype	Ba/F3 , 33%
------------------	-------------

Ba/F3 | 305224

Culture Medium RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3(820700a)

Supplements 5% FBS, 10ng/mL IL-3

Subculturing . 5 x 10⁵ /ml , 3 x 10⁵ ~ 1 x 10⁶ /ml .

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere 37°C, 5% CO₂ .

Flask Coating .

Freezing Procedure -78°C .

Ba/F3 | 305224

**Shipping
Conditions**

-78°C

**Storage
Conditions**

-150°C -196°C

.80°C

/ / HLA

Sterility

PCR