

hCMEC/D3 | 305024

Description	HCMEC/D3	, - (BBB)	.
	HCMEC/D3	, BBB	-5,
	HCMEC/D3	,	,

Organism	
Tissue	,

Synonyms	/D3, CMEC/D3,	/D3
-----------------	---------------	-----

Age	
------------	--

Gender	
---------------	--

Morphology	
-------------------	--

Cell type	
------------------	--

Growth properties	
--------------------------	--

Citation	hCMEC/D3(305024)
-----------------	-----------	---------

Biosafety level	1
------------------------	---

NCBI_TaxID	9606
-------------------	------

CellSaurusAccession	CVCL_U985
----------------------------	-----------

GMO Status	GMO-S1:	(hCMEC/D3)	SV40 T-	hTERT
-------------------	---------	------------	---------	-------

hCMEC/D3 | 305024

Viruses : 40(SV40)

Culture Medium EGM -2 MV -2 BulletKit(Lonza, Lonza CC-3202)

Supplements EBM-2

Freeze medium 50% + 40% FBS + 10% DMSO CM-1(Cytion 800100)

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere 37°C, 5% CO₂ .

Flask Coating

Freezing Procedure -78°C .

Shipping Conditions -78°C .

hCMEC/D3 | 305024

Storage Conditions	-150°C	-196°C	-80°C
---------------------------	--------	--------	-------

/ / HLA

Sterility	PCR
------------------	-----