

T2 | 305228

Description T2 T1 CTL T2 (TAP) T2

Organism

Synonyms T2(174 x CEM.T2), T2(174 x CEM.T2), 174xCEM.T2, CEMx721.174.T2

Morphology

Growth properties

Citation T2(305228)

Biosafety level 2

NCBI_TaxID 9606

CellosaurusAccession CVCL_2211

Culture Medium RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3(820700a)

Supplements 10% FBS

Subculturing : 22

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

T2 | 305228

Thawing and Culturing Cells									
1.									
2.		-150°C						3	
3.			37°C					40~60	
4.					70%				
5.			8ml			15ml			
6.		300 x g	3						
7.			10ml				T25		T25
8.									

Incubation Atmosphere

37°C, 5% CO₂

Flask Coating

Freezing Procedure

-78°C

Shipping Conditions

-78°C

Storage Conditions

-150°C -196°C .80°C

/ / HLA

Sterility

PCR