

Product sheet

NCI-H2170 | 305276

Description

NCI-H2170 . , . NCI-H2170
NCI-H2170 p63 . TP53 CDKN2A
NCI-H2170 .

Organism

Tissue

Disease

Synonyms

H2170, H-2170, NCH2170

Age

Gender

Ethnicity

Morphology

Growth
properties

Citation

NCI-H2170(305276)

Biosafety level

1

NCBI_TaxID

9606

CellosaurusAccession

CVCL_1535

Product sheet

NCI-H2170 | 305276

Culture Medium	RPMI 1640, w: 2.0mM	, w: 2.0g/L NaHCO3(	820700a)
----------------	---------------------	---------------------	----------

Supplements 10% FBS 2.5g/L 10mM HEPES

Dissociation Reagent	
----------------------	--

Subculturing PBS . T25 3~5ml, T75 5~10ml PBS . T25

Split ratio	1:3 1:6
-------------	---------

Fluid renewal 1~2

Freeze medium	(FBS) + 10% DMSO , CM-1(Cytion 800100)
---------------	---

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere	37°C, 5% CO2' .
-----------------------	-----------------

Flask Coating

Product sheet

NCI-H2170 | 305276

Freezing Procedure	-78°C	.	.
Shipping Conditions	-78°C	.	.
Storage Conditions	-150°C	-196°C	.80°C
Sterility	PCR	.	.