

Colo-320HSR | 305271

Description	COLO-320HSR	COLO-320
	COLO-320HSR 320HSR	(CEA)

Organism

Tissue

Disease

Synonyms	COLO320 HSR, COLO 320HSR, COLO 320 HSR
-----------------	--

Age	55
------------	----

Gender

Ethnicity

Morphology

Growth properties

Citation	COLO-320HSR(305271)
-----------------	----------------------

Biosafety level	1
------------------------	---

NCBI_TaxID	9606
-------------------	------

CellosaurusAccession	CVCL_1989
-----------------------------	-----------

Colo-320HSR | 305271

Protein expression , , , (ACTH),

Tumorigenic ,

Culture Medium RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3(820700a)

Supplements 10% FBS 2.5g/L 10mM HEPES

Dissociation Reagent

Subculturing PBS . T25 3~5ml, T75 5~10ml PBS . T25

Fluid renewal 2

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Colo-320HSR | 305271

Incubation Atmosphere 37°C, 5% CO₂

Flask Coating

Freezing Procedure -78°C

Shipping Conditions -78°C

Storage Conditions -150°C -196°C -80°C

/ / HLA

Sterility PCR