



**Bend.3 | 305265**

**Antigen expression** ICAM-1 +, VCAM-1 +, MAdCAM-1 +

**Viruses** : (A2 ) (MPyV) T (PyMT)

**Culture Medium** DMEM, w: 4.5g/L , w: 4mM L- , w: 3.7g/L NaHCO3, w: 1.0mM ( 820300a)

**Supplements** 10% FBS

**Dissociation Reagent**

**Subculturing** PBS . T25 3~5ml, T75 5~10ml PBS . T25

**Freeze medium** (FBS ) + 10% DMSO , CM-1(Cytion 800100)

**Thawing and Culturing Cells**

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

**Bend.3 | 305265**

**Incubation  
Atmosphere**

37°C, 5% CO<sub>2</sub>

**Flask Coating**

**Freezing  
Procedure**

-78°C

**Shipping  
Conditions**

-78°C

**Storage  
Conditions**

-150°C -196°C .80°C

/ / HLA

**Sterility**

PCR