

HEK293-CTLA4 | 305423

Description

HEK293-CTLA4 (APC) 4,000 ~ CTLA4 CD28 HEK293 inscreenex

Organism

Tissue

Disease / ; (HEK293)

Applications ;CTLA4 (,);T ; ;

Age

Gender

Morphology

Cell type

Growth properties

Citation HEK293-CTLA4(305423)

Biosafety level 1

NCBI_TaxID 9606

CellosaurusAccession CVCL_6G20

Product sheet

HEK293-CTLA4 | 305423

GMO Status	GMO-S1: HEK293	CTLA-4	.	.
Receptors expressed	CTLA4(CD152)			
Culture Medium	RPMI 1640, w: 2.0mM	, w: 2.0g/L NaHCO3(	820700a)	
Supplements	10% FBS, 1mM	, 10mM HEPES, 1% NEAA	. (G418-)	1mg/mL .
Dissociation Reagent	-EDTA			
Doubling time	24~36			
Subculturing	:	PBS	. PBS	CO2 /EDTA 2~3 T25
Split ratio	1~5			
Seeding density	$2\sim4 \times 10^4$ /cm ²			
Fluid renewal	2~3			
Post-Thaw Recovery	T25	1:2 ~ 1:3	24 ()	.
Freeze medium	(FBS) + 10% DMSO , CM-1(Cytion 800100)			

Product sheet

HEK293-CTLA4 | 305423

Thawing and
Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation
Atmosphere

37°C, 5% CO₂ ,

Flask Coating

Freezing
Procedure

-78°C .

Shipping
Conditions

-78°C .

Storage
Conditions

-150°C -196°C . 80°C .

Sterility

PCR .