

DI TNC1 | 305343

Description	
DI TNC1	1
DI TNC1	AMP(cAMP)
DI TNC1	B - (NF-κB)

Organism	
Tissue	

Disease	
Synonyms	
DITNC1, DI-TNC1, DI TNC-1	

Breed/Subspecies	

Age	
1	

Gender	

Morphology	
Cell type	
, II	

Growth properties	

Citation	
DI TNC1(305343)	

Biosafety level	
2	

NCBI_TaxID	
10116	

DI TNC1 | 305343

CellosaurusAccession CVCL_0247

GMO Status GMO-S1: (DI TNC1) GFAP SV40 .

Protein expression : 2 ,

Tumorigenic ,

Viruses : 40(SV40)

Culture Medium DMEM, w: 4.5g/L , w: 4mM L- , w: 3.7g/L NaHCO3, w: 1.0mM (820300a)

Supplements 10% FBS

Dissociation Reagent

Subculturing PBS . T25 3~5ml, T75 5~10ml PBS . T25

Fluid renewal 2~3

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

DI TNC1 | 305343

Thawing and Culturing Cells	
1.	.
2.	-150°C , 3 .
3.	37°C 40~60 .
4.	, 70% .
5.	8ml 15ml .
6.	300 x g 3 .
7.	10ml . T25 , T25
8.	. .

Incubation Atmosphere 37°C, 5% CO₂ .

Flask Coating

Freezing Procedure -78°C

Shipping Conditions -78°C

Storage Conditions -150°C -196°C .80°C

/ / HLA

Sterility PCR