

DOHH-2 | 305537

**Description**

DOHH-2 B B (DLBCL) .  
ABT-263 CHOP , .  
DOHH-2 . , CD20 Bcl-2 , MDM2

**Organism**

**Tissue**

**Disease**

B

**Synonyms**

DOHH2, DoHH-2, DOHH-2

**Age**

60

**Gender**

**Ethnicity**

**Growth  
properties**

**Citation**

DOHH-2 (Cytion 305537)

**Biosafety level**

1

**NCBI\_TaxID**

9606

**CellosaurusAccession**

CVCL\_1179

**GMO Status**

GMO-S1: (DOHH-2) EBV EBV , . B

Product sheet

DOHH-2 | 305537

|         |   |   |        |
|---------|---|---|--------|
| Viruses | : | - | (EBV). |
|---------|---|---|--------|

|                |                     |                     |          |
|----------------|---------------------|---------------------|----------|
| Culture Medium | RPMI 1640, w: 2.0mM | , w: 2.0g/L NaHCO3( | 820700a) |
|----------------|---------------------|---------------------|----------|

|             |     |     |
|-------------|-----|-----|
| Supplements | 10% | FBS |
|-------------|-----|-----|

|              |      |     |      |            |           |
|--------------|------|-----|------|------------|-----------|
| Subculturing | 15ml | PBS | (T25 | 3~5ml, T75 | 5~10ml ). |
|--------------|------|-----|------|------------|-----------|

|                 |                         |
|-----------------|-------------------------|
| Seeding density | 3-5×10 <sup>5</sup> /ml |
|-----------------|-------------------------|

|               |                   |   |             |         |
|---------------|-------------------|---|-------------|---------|
| Freeze medium | (FBS ) + 10% DMSO | , | CM-1(Cytion | 800100) |
|---------------|-------------------|---|-------------|---------|

|                             |    |         |   |       |     |     |   |
|-----------------------------|----|---------|---|-------|-----|-----|---|
| Thawing and Culturing Cells | 1. |         |   |       |     |     | . |
|                             | 2. | -150°C  |   |       | ,   | 3   | . |
|                             | 3. | 37°C    |   | 40~60 |     | .   |   |
|                             | 4. |         |   | ,     | 70% | .   |   |
|                             | 5. | 8ml     |   | 15ml  |     | .   |   |
|                             | 6. | 300 x g | 3 | .     |     |     |   |
|                             | 7. | 10ml    | . | T25   | ,   | T25 |   |
|                             | 8. | .       |   |       |     |     |   |

|                       |                           |   |
|-----------------------|---------------------------|---|
| Incubation Atmosphere | 37°C, 5% <sub>CO2</sub> , | . |
|-----------------------|---------------------------|---|

|                     |       |   |
|---------------------|-------|---|
| Shipping Conditions | -78°C | . |
|---------------------|-------|---|

Product sheet

DOHH-2 | 305537

Storage  
Conditions

-150°C -196°C .80°C

Sterility

PCR