

Product sheet

ND7/23 | 305520

Description	ND7/23 (DRG) (N18TG2) . , ,
	ND7/23 , (NGF) cAMP(db-cAMP) .
	ND7/23 TRPC , TRPC4 . TRPC4 RNA(shRNA)
Organism	,
Tissue	
Synonyms	ND7-23
Cell type	(N18 tg 2) ×
Growth properties	
Citation	ND7/23 (Cytion 305520)
Biosafety level	1
NCBI_TaxID	10090, 10116
CellosaurusAccession	CVCL_4259
Culture Medium	DMEM, w: 4.5g/L , w: 4mM L- , w: 3.7g/L NaHCO3, w: 1.0mM (820300a)
Supplements	10% FBS

Product sheet

ND7/23 | 305520

Seeding density	$1 - 3 \times 10^4 / \text{cm}^2$			
Freeze medium	(FBS) + 10% DMSO , CM-1(Cytion 800100)			

Thawing and Culturing Cells	1.			
	2.	-150°C	,	3
	3.	37°C	40~60	
	4.	,	70%	
	5.	8ml	15ml	
	6.	300 x g 3		
	7.	10ml	T25	T25
	8.			

Incubation Atmosphere	37°C, 5% _{CO2}
-----------------------	-------------------------

Shipping Conditions	-78°C
---------------------	-------

Storage Conditions	-150°C -196°C .80°C
--------------------	---------------------

Sterility	PCR
	,