

**Panc02-Luc | 305706**

**Description**  
Panc02-Luc Panc02 .Panc02 ,  
Panc02-Luc , , , .Panc02  
Panc02-Luc , , , , .

**Organism**  
**Tissue**

**Disease**  
**Synonyms** Panc02

**Breed/Subspecies** C57BL/6

**Age**  
**Gender**

**Growth properties**

**Citation** Panc02-Luc (Cytion 305706)

**Biosafety level** 1

**NCBI\_TaxID** 10090

**CellSaurusAccession** CVCL\_E31B

**Panc02-Luc | 305706**

<b>Protein expression</b>	Luc
---------------------------	-----

<b>Culture Medium</b>	RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3( 820700a)
-----------------------	--

<b>Supplements</b>	10% FBS
--------------------	---------

<b>Dissociation Reagent</b>	
-----------------------------	--

<b>Doubling time</b>	24~48
----------------------	-------

<b>Subculturing</b>	PBS . T25 3~5ml, T75 5~10ml PBS . T25
---------------------	---------------------------------------

<b>Seeding density</b>	$1\sim3 \times 10^4 / \text{cm}^2$
------------------------	------------------------------------

<b>Fluid renewal</b>	2~3
----------------------	-----

<b>Freeze medium</b>	+ 10% DMSO .
----------------------	--------------

<b>Thawing and Culturing Cells</b>	
1.	.
2.	-150°C , 3 .
3.	37°C 40~60 .
4.	, 70% .
5.	8ml 15ml .
6.	200 x g 5 .
7.	

Panc02-Luc | 305706

**Incubation Atmosphere** 37°C, 5% CO<sub>2</sub>

**Shipping Conditions** -78°C

**Storage Conditions** -150°C -196°C .80°C

/ / HLA