

HK/FDC | 300205

Description
HK/FDC HK
HK/FDC B , -IgM CD40 VCAM-1

Organism

Tissue

Disease

Applications B / B .FDC

Age

Gender

Ethnicity

Morphology

Cell type

Growth properties

Citation HK/FDC(300205)

Biosafety level 1

NCBI_TaxID 9606

HK/FDC | 300205

Viruses , i.A.

Culture Medium DMEM: F12(1:1), w: 3.1 g/L , w: 2.5mM L- , w: 15mM HEPES, w: 0.5mM , w: 1.2 g/L NaHCO3 (Cytion 820

Supplements 10% FBS

Dissociation Reagent

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere 37°C, 5% CO2' .

Flask Coating

Freezing Procedure -78°C .

HK/FDC | 300205

**Shipping
Conditions**

-78°C

**Storage
Conditions**

-150°C -196°C

.80°C

/ / HLA

Sterility

PCR