



**KU-19-19 | 305517**

**Mutational profile** : p.Glu17Lys,

**Culture Medium** RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3( 820700a)

**Supplements** 10% FBS

**Doubling time** ~48

**Freeze medium** + 10% DMSO

**Thawing and Culturing Cells**

- 1.
2. -150°C
3. 37°C 40~60
4. 70%
5. 8ml 15ml
6. 200 x g 5
- 7.

**Incubation Atmosphere** 37°C, 5% CO<sub>2</sub>

**Flask Coating**

**Shipping Conditions** -78°C

**Storage Conditions** -150°C -196°C .80°C

**KU-19-19 | 305517**

**/ / HLA**