

KU-19-19 | 305517

Mutational profile : p.Glu17Lys,

Culture Medium RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3(820700a)

Supplements 10% FBS

Doubling time ~48

Freeze medium + 10% DMSO

Thawing and Culturing Cells

- 1.
2. -150°C
3. 37°C 40~60
4. 70%
5. 8ml 15ml
6. 200 x g 5
- 7.

Incubation Atmosphere 37°C, 5% CO2

Flask Coating

Shipping Conditions -78°C

Storage Conditions -150°C -196°C .80°C

KU-19-19 | 305517

/ / HLA