

| 305850

Description (Pfeiffer) B (DLBCL) B ,
, B . BCL6 (NF-κB PI3K/AKT)
, Pfeiffer B , (ADCC), -CD20 , ,

Organism

Tissue

Disease B

Synonyms

Age

Gender

Ethnicity

Morphology

Cell type B

Growth properties

Citation (305850)

Biosafety level 1

NCBI_TaxID 9606

CellosaurusAccession CVCL_3326

| 305850

Mutational profile : p.Arg1171Cys,

Culture Medium RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3(820700a)

Supplements 10% FBS

Dissociation Reagent

Doubling time 24~30

Seeding density 2~10 × 10⁵/mL

Freeze medium + 10% DMSO .

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 200 x g 5 .
7. .

Incubation Atmosphere 37°C, 5% CO₂' .

| 305850

**Shipping
Conditions**

-78°C

**Storage
Conditions**

-150°C -196°C

.80°C

/ / HLA