

Sf9 | 604328

Description	Sf9	Sf21 (IPLB-Sf-21-AE)	.
	Sf9		. Sf9
	1	26~28	. Sf9 / Sf9
	Sf9	Sf21	, Sf21 Sf9 .
	Sf9	Spodoptera frugiperda rhabdovirus(SfRV)	, Sf9 . Sf9

Organism

Tissue

Applications

Synonyms Sf9, sf9, SF-9, Sf-9, sf-9, Sf 9, Spodoptera frugiperda 9, Sf 9, IPLB-Sf-9AE, IPLB-SF-9AE, IPLB-SF-9, IPLB-Sf-9, IPLB-Sf9, IPLB-Sf9

Age

Gender

Morphology

Growth properties

Citation Sf9(604328)

Biosafety level 1

NCBI_TaxID 7108

CellosaurusAccession CVCL_0549

Sf9 | 604328

Virus susceptibility , (MNPV), (SLE)

Culture Medium (PAN)

Supplements 2% FBS

Dissociation Reagent

Subculturing . 15ml . 5ml

Seeding density 1×10^4 /cm². 26~30 , .

Fluid renewal 2~3

Freeze medium (FBS) + 10% DMSO , CM-1(800100)

Sf9 | 604328

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere

27°C, 0% CO₂, humidified atmosphere.

Shipping Conditions

Cryopreserved cell lines are shipped on dry ice in validated, insulated packaging with sufficient refrigerant to maintain approximately -78 °C throughout transit. On receipt, inspect the container immediately and transfer vials without delay to appropriate storage.

Storage Conditions

For long-term preservation, place vials in vapor-phase liquid nitrogen at about -150 to -196 °C. Storage at -80 °C is acceptable only as a short interim step before transfer to liquid nitrogen.

/ / HLA

Sterility

PCR .