

CHO-PDCD1 | 305973

Description	: . 6,250 . - . , CHO-PDCD1 T , B 1(PD-1; PDCD1/CD279) PDCD1 CHO . CHO-PDCD1 , , - , .		
Organism			
Tissue			
Disease	, ; PD-1 (PDCD1/CD279)		
Applications	; PD-1 ; ; T ;		
Age			
Gender			
Morphology			
Cell type			
Growth properties	/		
Citation	CHO-PDCD1 (Cytion 305973)		
Biosafety level	1		
NCBI_TaxID	10029		

Product sheet

CHO-PDCD1 | 305973

CellosaurusAccession CVCL_A8X0

GMO Status GMO-S1: CHO PDCD1 . , .

Receptors expressed PDCD1/CD279

Culture Medium : DMEM: F12(1:1), w: 3.1g/L , w: 2.5mM L- , w: 15mM HEPES, w: 0.5mM , w: 1.2g/L NaHCO3(Cytion
CHO A(InSCREENeX, InSCREENeX INS-ME-1039)

Supplements : 5% FBS . (G418-) 0.5mg/mL .

Dissociation Reagent : -EDTA

Doubling time 14~16

Subculturing : PBS . PBS CO2 /EDTA (: T25~3

Split ratio 1~5

Seeding density $2\sim5 \times 10^4$ /cm²

Fluid renewal 2~3

Post-Thaw Recovery T25 1:2 ~ 1:3 24 () .

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

CHO-PDCD1 | 305973

Thawing and
Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation
Atmosphere

37°C, 5% CO₂ ,

Shipping
Conditions

-78°C .

Storage
Conditions

-150°C -196°C . 80°C .

Sterility

PCR .
, .