

Description		:	/	.	6,250	.
	-	*	,			
CHO-CD36	CD36	(CHO)	. CD36	IV(GPIV)	(FAT)	
CHO-CD36	-	, ,	CD36	.	,	, ,

Organism

Disease	CD36

Applications	CD36			

Gender

Morphology

[illegible]

Growth properties

Citation	CHO-CD36 (Cytion	305979)
-----------------	------------------	---------

Biosafety level	1
------------------------	---

NCBI_TaxID	10029
------------	-------

Product sheet

CHO-CD36 | 305979

CellosaurusAccession CVCL_8848

GMO Status GMO-S1: CHO CD36 . , .

Receptors expressed CD36

Culture Medium : DMEM: F12(1:1), w: 3.1g/L , w: 2.5mM L- , w: 15mM HEPES, w: 0.5mM , w: 1.2g/L NaHCO3(Cytion
CHO A(InSCREENeX, InSCREENeX INS-ME-1039)

Supplements : 5% FBS . (G418-) 0.5mg/mL .

Dissociation Reagent : -EDTA

Doubling time 14~16

Subculturing : PBS . PBS CO2 /EDTA (: T25~3

Split ratio 1~5

Seeding density $2\sim5 \times 10^4$ /cm²

Fluid renewal 2~3

Post-Thaw Recovery T25 1:2 ~ 1:3 24 () .

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

CHO-CD36 | 305979

Thawing and
Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation
Atmosphere

37°C, 5%_{CO2} .

Shipping
Conditions

-78°C .

Storage
Conditions

-150°C -196°C . 80°C .

Sterility

PCR .
, .