

4T1-Luc | 305663

Description	4T1-Luc	,	4T1	.	4T1	,4
		(BLI)	,	.	.	
	4T1-Luc			.	,	,

Organism
Tissue

Disease

Breed/Subspecies	BALB/cfC3H
-------------------------	------------

Gender

Morphology

Growth properties

Citation	4T1-Luc (Cytion 305663)
-----------------	-------------------------

Biosafety level	1
------------------------	---

NCBI_TaxID	10090
-------------------	-------

CellosaurusAccession	CVCL_J239
-----------------------------	-----------

Antigen expression	Luc
---------------------------	-----

4T1-Luc | 305663

Tumorigenic , BALB/c .

MSI-status

Culture Medium RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3(820700a)

Supplements 10% FBS

Dissociation Reagent

Subculturing PBS . T25 3~5ml, T75 5~10ml PBS . T25

Seeding density $1\sim3 \times 10^4 / \text{cm}^2$

Fluid renewal 2~3

Freeze medium + 10% DMSO .

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 200 x g 5 .
7. .

4T1-Luc | 305663

**Incubation
Atmosphere**

37°C, 5% CO₂

**Shipping
Conditions**

-78°C

**Storage
Conditions**

-150°C -196°C .80°C

/ / HLA