

DAKIKI | 305928

Description

DAKIKI IgA1 IgA1 - (EBV) B-
NeuAc:GalNAc-IgA1 α 2,6- , DAKIKI IgA1 O-
, DAKIKI B B . DAKIKI IgA

Organism

Tissue

Synonyms , 1

Age

Gender

Ethnicity

Morphology

Cell type B

Growth properties

Citation DAKIKI (Cytion 305928)

Biosafety level 2

NCBI_TaxID 9606

CellosaurusAccession CVCL_3675

Product sheet

DAKIKI | 305928

Culture Medium	RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3(820700a)
Supplements	10% FBS
Dissociation Reagent	
Subculturing	⁵ /ml , , ml .. 1 x 10
Seeding density	2~5 × 10 ⁵ /cm ²
Fluid renewal	2~3
Freeze medium	+ 10% DMSO .
Thawing and Culturing Cells	1. . 2. -150°C , 3 . 3. 37°C 40~60 . 4. , 70% . 5. 8ml 15ml . 6. 200 x g 5 . 7.
Incubation Atmosphere	37°C, 5% _{CO2} , .

Product sheet

DAKIKI | 305928

Shipping
Conditions

-78°C

Storage
Conditions

-150°C

-196°C

.80°C