

HEK293-Rpn11-HTBH | 305719

Description

HEK293 -Rpn11-HTBH HEK293 (Human Embryonic Kidney 293) ,26S
JAMM , . HTBH (-TEV- -)
26S , - (UPS) , Rpn11/PSMD14 , ,

Organism

Tissue

Disease / (HEK293 ; Rpn11-HTBH)

Applications 26S ; Rpn11/PSMD14 ; - ; ; ; HTBH- ;

Morphology

Cell type

Growth properties

Citation HEK293 -Rpn11-HTBH (Cytion 305719)

Biosafety level 1

NCBI_TaxID 9606

GMO Status GMO-S1: HEK293 Rpn11-HTBH (-TEV- - Rpn11/PSMD14)

Culture Medium DMEM, w: 4.5g/L , w: 4mM L- , w: 3.7g/L NaHCO3, w: 1.0mM (820300a)

Product sheet

HEK293-Rpn11-HTBH | 305719

Supplements 10% FBS

Dissociation Reagent PBS , 5 Accutase .

Doubling time 24~36

Subculturing , PBS , Accutase 8~10 . 300×g 3 ,

Split ratio 1 10

Seeding density $2\sim4 \times 10^4$ /cm²

Fluid renewal 2~3

Freeze medium + 10% DMSO .

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 200 x g 5 .
7. .

Incubation Atmosphere 37°C, 5% CO₂ ,

Shipping Conditions -78°C .

Product sheet

HEK293-Rpn11-HTBH | 305719

Storage
Conditions

-150°C -196°C

.80°C

.