

Product sheet

Hs 611. T | 305915

Description

Hs 611.T 48 (cHL) . cHL T
Hs 611.T HRS , NF-κB , JAK-STAT , (CD30, CD15, CD25) (anti-CD30 ADC)

Organism

Tissue

Disease

Metastatic site

Applications ; NF-κB JAK-STAT ; CD30/CD15/CD25 ; ;

Synonyms Hs-611-T, Hs611T, HS611T

Age 48

Gender

Ethnicity

Morphology

Cell type B

Growth properties

Citation Hs 611.T (Cytion 305915)

Biosafety level 1

NCBI_TaxID 9606

Product sheet

Hs 611. T | 305915

CellosaurusAccession CVCL_0823

Culture Medium DMEM, w: 4.5g/L , w: 4mM L- , w: 3.7g/L NaHCO3, w: 1.0mM (820300a)

Supplements 10% FBS

Dissociation Reagent

Doubling time 24~48

Subculturing $3-5 \times 10^5$ /ml . (300×g, 5) .

Split ratio 1~3

Seeding density $2-4 \times 10^5$ /ml

Fluid renewal 2~3

Freeze medium + 10% DMSO .

Product sheet

Hs 611. T | 305915

Thawing and
Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 200 x g 5 .
7. .

Incubation
Atmosphere

37°C, 5%_{CO2} .

Shipping
Conditions

-78°C .

Storage
Conditions

-150°C -196°C . 80°C .