

HROC60 | 300827

**Description** 2006 PD CRC .

**Organism**

**Tissue** , UICC I

**Disease** , TNM T2N0M0R0L0V0, G2, Lk(n) + 0,  $\Sigma$  Lk(n) 47

**Age** 71

**Gender**

**Ethnicity**

**Morphology**

**Growth properties**

**Citation** HROC60( 300827)

**Biosafety level** 1

**NCBI\_TaxID** 9606

**CellSaurusAccession** CVCL\_1G05

**Protein expression** PTEN

**Antigen expression** Her2/neu +, EGFR +, CD326+, CD44+, CD54+, CD47+, CD71+, CD15-, CD73+, CD95+, CD274+, CD133low , CD276+, IDOweak, MHC-I+, CD133weak, CD66acdewweak, EpCAM+, IFN- $\gamma$  MHCII+ , cFLIPweak

**HROC60 | 300827**

**Tumorigenic** ,

**Viruses** SV40, JC/BK, HBV, HCV, HIV .

**Ploidy status**

**MSI-status** MSS

**Mutational profile** KRAS

**Culture Medium** DMEM: F12(1:1), w: 3.1 g/L , w: 2.5mM L- , w: 15mM HEPES, w: 0.5mM , w: 1.2 g/L NaHCO3 (Cytion 820

**Supplements** 10% FBS

**Dissociation Reagent**

**Doubling time** 29

**Subculturing** PBS . T25 3~5ml, T75 5~10ml PBS . T25

**Seeding density**  $2 \times 10^4$  /cm<sup>2</sup>

**Fluid renewal** 3~5

**Post-Thaw Recovery** 1~2

**Freeze medium** (FBS ) + 10% DMSO , CM-1(Cytion 800100)

HROC60 | 300827

**Thawing and Culturing Cells**

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

**Incubation Atmosphere**

37°C, 5% CO<sub>2</sub> .

**Flask Coating**

**Freezing Procedure**

-78°C .

**Shipping Conditions**

-78°C .

**Storage Conditions**

-150°C -196°C .80°C .

/ / HLA

**Sterility**

PCR .