

HK-CRISPR-Nup93-mEGFP | 300655

Description	HK-CRISPR-Nup93-mEGFP	, CRISPR/Cas9	(mEGFP)	Nup93
	Nup93	. MEGFP	.	.
	HK-CRISPR-Nup93-mEGFP	,	.	.

Organism
Tissue

Disease

Age	30
------------	----

Gender
Ethnicity

Morphology
Growth properties

Citation	HK-CRISPR-Nup93-mEGFP(300655)
-----------------	------------------------	---------

Biosafety level	1
NCBI_TaxID	9606

Depositor	(EMBL)				
GMO Status	GMO-S1: HeLa Kyoto	Nup93	mEGFP	(knock-in)	.

HK-CRISPR-Nup93-mEGFP | 300655

Protein expression Nup153, mEGFP

Culture Medium DMEM, w: 4.5g/L , w: 4mM L- , w: 3.7g/L NaHCO3, w: 1.0mM (820300a)

Supplements 10% FBS

Dissociation Reagent

Subculturing PBS . T25 3~5ml, T75 5~10ml PBS . T25

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

- Thawing and Culturing Cells**
1. .
 2. -150°C , 3 .
 3. 37°C 40~60 .
 4. , 70% .
 5. 8ml 15ml .
 6. 300 x g 3 .
 7. 10ml . T25 , T25
 8. .

Incubation Atmosphere 37°C, 5% CO2' .

Flask Coating

HK-CRISPR-Nup93-mEGFP | 300655

Freezing Procedure	-78°C	.
---------------------------	-------	---

Shipping Conditions	-78°C	.
----------------------------	-------	---

Storage Conditions	-150°C	-196°C	.80°C	.
---------------------------	--------	--------	-------	---

/ / HLA

Sterility	PCR	.
------------------	-----	---