

L-540 | 300201

Description	L-540	.	B						
	L-540			CD30	.	L-540	-		CD30

Organism

Tissue

Disease

Synonyms	L 540, L540
-----------------	-------------

Age	20
------------	----

Gender

Ethnicity

Morphology

Growth properties

Citation	L-540(300201)
-----------------	--------	---------

Biosafety level	1
------------------------	---

NCBI_TaxID	9606
-------------------	------

CellosaurusAccession	CVCL_1362
-----------------------------	-----------

L-540 | 300201

Viruses	EBV
----------------	-----

Culture Medium	RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3(820700a)
-----------------------	--

Supplements	10% FBS
--------------------	---------

Subculturing	⁵ /ml , , ml .. 1 x 10
---------------------	-----------------------------------

Fluid renewal	3
----------------------	---

Freeze medium	(FBS) + 10% DMSO , CM-1(Cytion 800100)
----------------------	---

Thawing and Culturing Cells	1. .
	2. -150°C , 3 .
	3. 37°C 40~60 .
	4. , 70% .
	5. 8ml 15ml .
	6. 300 x g 3 .
	7. 10ml . T25 , T25
	8. .

Incubation Atmosphere	37°C, 5% CO ₂ ' .
------------------------------	------------------------------

Flask Coating	.
----------------------	---

L-540 | 300201

Freezing Procedure	-78°C	.
---------------------------	-------	---

Shipping Conditions	-78°C	.
----------------------------	-------	---

Storage Conditions	-150°C	-196°C	.80°C	.
---------------------------	--------	--------	-------	---

/ / HLA

Sterility	PCR	.
------------------	-----	---