

Product sheet

CLS-138 | 400177

Description

	NMRI	CLS-138	.	,	
CLS-138			.		
CLS-138			.		
NMRI	CLS-138			.	,
CLS-138				.	NMRI

Organism

Tissue

Disease

Breed/Subspecies NMRI

Age

Gender

Morphology

Cell type

Growth
properties

Citation CLS-138(400177)

Biosafety level 1

NCBI_TaxID 10090

CellosaurusAccession CVCL_5726

CLS-138 | 400177

Tumorigenic ,

Culture Medium DMEM, w: 4.5g/L , w: 4mM L- , w: 3.7g/L NaHCO₃, w: 1.0mM (820300a)

Supplements 10% FBS

Dissociation Reagent

Subculturing PBS . T25 3~5ml, T75 5~10ml PBS . T25

Seeding density 2×10^4 /cm² 2 .

Fluid renewal 3~5

Post-Thaw Recovery 5×10^4 , 24 .

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

Product sheet

CLS-138 | 400177

Thawing and
Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation
Atmosphere

37°C, 5% CO₂ ,

Flask Coating

Freezing
Procedure

-78°C .

Shipping
Conditions

-78°C .

Storage
Conditions

-150°C -196°C . 80°C .

Sterility

PCR .