

Yac-1 | 400383

Description	(Mo-MuLV)	.	(NK)	NK	.
Organism					
Tissue					
Disease					
Synonyms	YAC-1, YAC				
Breed/Subspecies	A/Sn				
Cell type					
Growth properties					
Citation	Yac-1(400383)				
Biosafety level	1				
NCBI_TaxID	10090				
CellosaurusAccession	CVCL_2244				
Culture Medium	RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3(820700a)				
Supplements	10% FBS				

Yac-1 | 400383

Subculturing		5×10^5 /ml	,	$3 \times 10^5 \sim 1 \times 10^6$ /ml	
Fluid renewal	2~3				
Freeze medium		(FBS) + 10% DMSO	,	CM-1(Cytion 800100)	
Thawing and Culturing Cells	1.				
	2.	-150°C	,	3	
	3.	37°C		40~60	
	4.		,	70%	
	5.	8ml		15ml	
	6.	300 x g	3		
	7.	10ml		T25	, T25
	8.				
Incubation Atmosphere		37°C, 5% CO ₂			
Flask Coating					
Freezing Procedure		-78°C			
Shipping Conditions		-78°C			
Storage Conditions		-150°C -196°C		.80°C	

/ / HLA

Yac-1 | 400383

Sterility

PCR