

HT22 | 305158

| | | | | | | |
|--------------------|------|-----|---|---|--------|------|
| Description | HT22 | HT4 | , | . | SV40 T | HT22 |
| | HT22 | , | , | . | | |
| | HT22 | | | | | |
| | HT22 | , | | . | | |

Organism

| | |
|---------------|---|
| Tissue | , |
|---------------|---|

Synonyms HT-22

Morphology

Growth properties

| | | |
|-----------------|-------|---------|
| Citation | HT22(| 305158) |
|-----------------|-------|---------|

Biosafety level 1

| | |
|-------------------|-------|
| NCBI_TaxID | 10090 |
|-------------------|-------|

CellosaurusAccession CVCL_0321

| | | | | |
|-------------------|---------|--------|--------|---|
| GMO Status | GMO-S1: | (HT22) | SV40 T | . |
|-------------------|---------|--------|--------|---|

| | | | | | |
|-----------------------|-----------------|-------------|---|---|----------|
| Culture Medium | DMEM, w: 4.5g/L | , w: 4mM L- | , w: 3.7g/L NaHCO ₃ , w: 1.0mM | (| 820300a) |
|-----------------------|-----------------|-------------|---|---|----------|

HT22 | 305158

Supplements 10% FBS

Dissociation Reagent

Subculturing PBS . T25 3~5ml, T75 5~10ml PBS . T25

Freeze medium 50% + 40% FBS + 10% DMSO CM-1(Cytion 800100)

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere 37°C, 5% CO₂ .

Flask Coating

Freezing Procedure -78°C .

Shipping Conditions -78°C .

Product sheet

HT22 | 305158

| | | | |
|---------------------------|--------|--------|-------|
| Storage Conditions | -150°C | -196°C | -80°C |
|---------------------------|--------|--------|-------|

/ / HLA

| | |
|------------------|-----|
| Sterility | PCR |
|------------------|-----|