

MR1 | 305000

Description	MR1	NS-1	,	T	,	Th1	.	,	CD40	
	MR1	CD154	CD40	,					-CD154	T
	MR1	CD40-CD154	T	B					. MR1	

Organism

Morphology

Growth properties

Citation	MR1(305000)
-----------------	------	---------

Biosafety level 1

NCBI_TaxID	10090/10032
-------------------	-------------

CellosaurusAccession CVCL_8964

Protein expression	,	,	CD40	(CD154, CD40L, gp39)
---------------------------	---	---	------	----------------------

Culture Medium	RPMI 1640, w: 2.0mM	,	w: 2.0g/L NaHCO3(820700a)
-----------------------	---------------------	---	-------------------	----------

Supplements 10% FBS, 0.05mM 2-

MR1 | 305000

Subculturing 5 /ml , , ml .. 1 x 10

Fluid renewal 2~3

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere 37°C, 5% CO₂ .

Flask Coating

Freezing Procedure -78°C .

Shipping Conditions

-78°C .

Storage Conditions -150°C -196°C . 80°C .

/ / HLA

Product sheet

MR1 | 305000

Sterility	PCR
------------------	-----