

MR1 | 305000

Description	MR1	NS-1	,	T	,	Th1	.	,	CD40	
	MR1	CD154	CD40	,					-CD154	T
	MR1	CD40-CD154	T	B					. MR1	

Organism

Morphology

Growth properties

Citation	MR1(305000)
-----------------	------	---------

Biosafety level 1

NCBI_TaxID	10090/10032
-------------------	-------------

CellosaurusAccession CVCL_8964

Protein expression	,	,	CD40	(CD154, CD40L, gp39)
---------------------------	---	---	------	----------------------

Culture Medium	RPMI 1640, w: 2.0mM	,	w: 2.0g/L NaHCO3(820700a)
-----------------------	---------------------	---	-------------------	----------

Supplements 10% FBS, 0.05mM 2-

MR1 | 305000

Subculturing	5	/ml	,	,	ml	..	1 x 10
Fluid renewal	2~3						
Freeze medium				(FBS) + 10% DMSO	,		CM-1(Cytion 800100)
Thawing and Culturing Cells	1.						
	2.	-150°C	,		3		
	3.			37°C	40~60		
	4.		,	70%			
	5.			8ml	15ml		
	6.	300 x g	3				
	7.	10ml			T25	,	T25
	8.						
Incubation Atmosphere	37°C, 5%	CO ₂					
Flask Coating							
Freezing Procedure		-78°C					
Shipping Conditions		-78°C					
Storage Conditions		-150°C	-196°C		.80°C		

Product sheet

MR1 | 305000

/ / HLA

Sterility	PCR	.
------------------	-----	---