

B-LCL-HROC39(Bc HROC39) | 300863

Description
B-LCL-HROC39 (Bc HROC39) B - (EBV) B .
, B-LCL-HROC39(Bc HROC39) CD19 CD20 CD23, CD80 B
, B-LCL-HROC39(Bc HROC39) (: IgG, IgM IgA) , .

Organism
Tissue

Disease
Synonyms B-LCL CO39, Bc HROC39

Age 69

Gender
Ethnicity

Morphology
Cell type B

Growth properties

Citation B-LCL-HROC39(Bc HROC39)(300863)

Biosafety level 2

NCBI_TaxID 9606

B-LCL-HROC39(Bc HROC39) | 300863

CellosaurusAccession CVCL_YD56

Surface antigens CD19

Viruses : EBV

Culture Medium RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3(820700a)

Supplements 10% FBS

Subculturing ⁵ /ml , , ml .. 1 x 10

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere 37°C, 5% CO2' .

B-LCL-HROC39(Bc HROC39) | 300863

Flask Coating

Freezing Procedure -78°C

Shipping Conditions -78°C

Storage Conditions -150°C -196°C .80°C

/ / HLA

Sterility PCR