

- (HMSC-BM) | 300665

Description	(HMSC-BM)				
	MSC		HMSC-BM		
	HMSC-BM	6	92%~95%	1 x 10	

Organism

Tissue

Applications

Age

Gender

Ethnicity

Morphology

Cell type

Growth
properties

Citation

Biosafety level

NCBI_TaxID

Product sheet

- (HMSC-BM) | 300665

Antigen expression	CD73/CD90/CD105() CD14/CD34/CD45/HLA-DR()		MSC(P2-P3)	
Viruses	HBV(PCR),	(PCR), HIV-1/2(IFA)	HBV, HCV, HSV1, HSV2, CMV, EBV, HHV6,	
Culture Medium	MEM, w: 2.0mM	, :	, w: 1.0mM	, w: 2.2g/L NaHCO3
Supplements	10% FBS, 2ng/mL bFGF			
Dissociation Reagent	-EDTA			
Subculturing	:	PBS	. PBS	CO2 /EDTA 2-3 T25
Seeding density	$1\sim3 \times 10^4 /cm^2$			
Fluid renewal	24	2~3		
Freeze medium	80% FBS + 10%	+ 10% DMSO	,	CM-1(800100)

- (HMSC-BM) | 300665

Thawing and
Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation
Atmosphere

37°C, 5%_{CO2} .

Flask Coating

Freezing
Procedure

-78°C .

Shipping
Conditions

-78°C .

Storage
Conditions

-150°C -196°C . 80°C .

Sterility

PCR .