

MH-TC-5123 | 500287

**Description**

**Organism**

**Tissue**

**Disease**

**Breed/Subspecies**

**Age** 10

**Gender**

**Morphology**

**Growth properties**

**Citation** MH-TC-5123( 500287)

**Biosafety level** 1

**NCBI\_TaxID** 10116

**CellosaurusAccession** CVCL\_5800

**Tumorigenic** ,

**Viruses** RAP .

MH-TC-5123 | 500287

**Culture Medium** RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3( 820700a)

**Supplements** 10% FBS

**Dissociation Reagent**

**Subculturing** PBS . T25 3~5ml, T75 5~10ml PBS . T25

**Seeding density**  $1 \times 10^4 / \text{cm}^2$

**Fluid renewal** 3~5

**Freeze medium** (FBS ) + 10% DMSO , CM-1(Cytion 800100)

- Thawing and Culturing Cells**
1. .
  2. -150°C , 3 .
  3. 37°C 40~60 .
  4. , 70% .
  5. 8ml 15ml .
  6. 300 x g 3 .
  7. 10ml . T25 , T25
  8. .

**Incubation Atmosphere** 37°C, 5% CO<sub>2</sub> .

**Flask Coating**

MH-TC-5123 | 500287

<b>Freezing Procedure</b>	-78°C	.
---------------------------	-------	---

<b>Shipping Conditions</b>	-78°C	.
----------------------------	-------	---

<b>Storage Conditions</b>	-150°C	-196°C	.80°C	.
---------------------------	--------	--------	-------	---

/ / HLA

<b>Sterility</b>	PCR	.
------------------	-----	---