

DLD-1 | 300220

<b>Description</b>	DLD-1	.	.	DLD-1
		13	KRAS	. DLD-1 APC

**Organism**

**Tissue**

**Disease**

<b>Synonyms</b>	DLD 1, DLD1, CoCL3
-----------------	--------------------

<b>Age</b>	67
------------	----

**Gender**

**Morphology**

**Growth properties**

<b>Citation</b>	DLD-1(	300220)
-----------------	--------	---------

<b>Biosafety level</b>	1
------------------------	---

<b>NCBI_TaxID</b>	9606
-------------------	------

<b>CellosaurusAccession</b>	CVCL_0248
-----------------------------	-----------

**Protein expression**

DLD-1 | 300220

**Tumorigenic**

**Viruses**

**Products** (CEA) 0.5ng/10 exp6 /10 ,

**Karyotype** 2n = 46

**Culture Medium** RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3( 820700a)

**Supplements** 10% FBS

**Dissociation Reagent**

**Doubling time** 15

**Subculturing** PBS . T25 3~5ml, T75 5~10ml PBS . T25

**Seeding density**  $1\sim 2 \times 10^4 / \text{cm}^2$

**Fluid renewal** 2~3

**Freeze medium** (FBS ) + 10% DMSO , CM-1(Cytion 800100)

DLD-1 | 300220

Thawing and Culturing Cells	
1.	.
2.	-150°C , 3 .
3.	37°C 40~60 .
4.	, 70% .
5.	8ml 15ml .
6.	300 x g 3 .
7.	10ml . T25 , T25
8.	. .

**Incubation Atmosphere** 37°C, 5% CO<sub>2</sub> .

**Flask Coating**

**Freezing Procedure** -78°C . . . . .

**Shipping Conditions** -78°C . . . . .

**Storage Conditions** -150°C -196°C .80°C . . . . .

/ / HLA

**Sterility** PCR . . . . .