

RCC-LR | 300236

**Description** 1999 63 pT3a, No, M1/GIII . HLA-A2.1 .

**Organism**

**Tissue**

**Disease**

**Synonyms** KTCTL120, RCC-LR

**Age** 63

**Gender**

**Ethnicity**

**Morphology**

**Growth properties** ,

**Citation** RCC-LR( 300236)

**Biosafety level** 1

**NCBI\_TaxID** 9606

**CellosaurusAccession** CVCL\_5883

**Surface antigens** 8,18,19,

RCC-LR | 300236

<b>Protein expression</b>	IL8
<b>Mutational profile</b>	IL8 RS1126647 3-UTR Wt
<b>Culture Medium</b>	RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3( 820700a)
<b>Supplements</b>	10% FBS
<b>Dissociation Reagent</b>	
<b>Subculturing</b>	PBS . T25 3~5ml, T75 5~10ml PBS . T25
<b>Split ratio</b>	1:2 ~ 1:3
<b>Fluid renewal</b>	1~2
<b>Freeze medium</b>	(FBS ) + 10% DMSO , CM-1(Cytion 800100)

RCC-LR | 300236

**Thawing and Culturing Cells**

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

**Incubation Atmosphere**

37°C, 5% CO<sub>2</sub> .

**Flask Coating**

**Freezing Procedure**

-78°C .

**Shipping Conditions**

-78°C .

**Storage Conditions**

-150°C -196°C .80°C .

/ / HLA

**Sterility**

PCR .

RCC-LR | 300236

---

**STR**

**Amelogenin:** x,x  
**CSF1PO:** 12  
**D13S317:** 12, 14  
**D16S539:** 12  
**D5S818:** 13  
**D7S820:** 11, 12  
**TH01:** 7, 8  
**TPOX:** 8.1  
**vWA:** 16, 17  
**D3S1358:** 16, 17  
**D21S11:** 29.3  
**D18S51:** 13, 14  
**Penta E:** 12  
**Penta D:** 9,14  
**D8S1179:** 14, 15  
**FGA:** 20,22  
**PEZ6:** HROG15