

MMQ | 300498

Viruses	SMRV-		
Products			
Karyotype	6% - 49-522n> -		
Culture Medium	RPMI 1640, w: 2.0mM	, w: 2.0g/L NaHCO3(820700a)
Supplements	7.5%	, 2.5%	FBS
Subculturing		5 x 10 ⁵ /ml	, 3 x 10 ⁵ ~ 1 x 10 ⁶ /ml
Seeding density	2 x 10 ⁵ /mL		
Freeze medium		(FBS) + 10% DMSO	, CM-1(Cytion 800100)
Thawing and Culturing Cells	<ol style="list-style-type: none"> -150°C 37°C 70% 8ml 15ml 300 x g 3 10ml T25 		
Incubation Atmosphere	37°C, 5% CO ₂		

MMQ | 300498

Flask Coating

Freezing Procedure

-78°C

Shipping Conditions

-78°C

Storage Conditions

-150°C -196°C .80°C

/ / HLA

Sterility

PCR